

AGP

30 days | Tue Jun 10, 2025 - Wed Jul 9, 2025

180

Time in Ranges Goals for Type 1 and Type 2 Diabetes

Each 5% increase in the Target Range is clinically beneficial.
Each 1% time in range = about 15 minutes per day




Target Range: 70-180 mg/dL
Very High: Above 250 mg/dL
Very Low: Below 54 mg/dL

Glucose Metrics

Average Glucose Goal: <154 mg/dL	135 mg/dL
GMI Goal: <7%	6.5%
Coefficient of Variation Goal: <36%	20.0%
Time CGM Active	98.6%


Ambulatory Glucose Profile (AGP)

AGP is a summary of glucose values from the report period, with median (50%) and other percentiles shown as if they occurred in a single day.



Daily Glucose Profile

Each daily profile represents a midnight-to-midnight period.



Patent pending - HealthPartners Institute dba International Diabetes Center - All Rights Reserved. ©2022